Minnesota Department of Health Radiation Control, X-ray Unit

Processor Quality Control/Step Wedge Test

Room #											
Step Wedge Technique mA					Time		kVp1	Distance			
Date	Temperature	-1	0	+1	Initial	Date	Temperature	-1	0	+1	Initial
L	I .	1				l		1	ı	1	

If Step Wedge Test is OFF by two steps (+ or -), corrective action must be taken and stepwedge test repeated if processor is in control. Standard step_____.



Radiation Control, X-ray Unit 625 North Robert Street PO Box 64497 St. Paul, Minnesota 55164-0497 651-201-4545 www.health.state.mn.us/xray health.xray@state.mn.us