Minnesota Department of Health Radiation Control, X-ray Unit

Repeat/Reject Analysis Worksheet

Time Period:	Time Period: to		
Total # Images Used:			
Cause	Number of Images	Total Number	Total Percentage
Image - Black			
Image - Dark			
Image- Good			
Image - Light			
Fog - Cassette			
Fog - Darkroom			
Mechanical			
Other			
Patient Motion			
Positioning			
Static			
-	lysis must include at a minimum the nary of causes for the retakes. Inclu		•
<u> </u>	Repeat Rate = Total Number Repeats Total Images Taken Quarterly Repeat Rate		
Optional: Individu	ıal Cause Repeat Rate = ' <u>Total Num</u> b	per Repeats for That (Cause
1	Total Number of Repeats x 100		



Radiation Control, X-ray Unit 625 North Robert Street PO Box 64497 St. Paul, Minnesota 55164-0497 651-201-4545 www.health.state.mn.us/xray