

Pure Tone Audiometry Worksheet “Sweep Screening”

FOR CHILDREN 3 YEARS THROUGH 10 YEARS

Conditioning

Place headphones on the table facing the child. Present 2000 Hz at maximum volume. Practice raising either hand when a tone is heard.

Place headphones on the child, red headphone on the child’s right ear and blue headphone on the child’s left ear. Present a conditioning tone of 40dB at 1000 Hz.

Sweep Screening

Right: 20 dB at 1000, 2000, 4000 Hz

Switch to

Left: 20 dB at 4000, 2000, 1000 Hz

Increase to 25dB at 500 Hz

Switch to **Right:** 25 dB at 500 Hz

Present each sound 1-2 times and vary the spacing or timing between sounds.

Volume (dB)	25	20	20	20
Frequency (Hz)	500	1000	2000	4000
Right Ear	Finish	Start	→	
Left Ear	↑	←		↓

Pure Tone Audiometry Worksheet “Sweep Screening” FOR CHILDREN 3 YEARS THROUGH 10 YEARS

Conditioning

Place headphones on the table facing the child. Present 2000 Hz at maximum volume. Practice raising either hand when a tone is heard.

Place headphones on the child, red headphone on the child’s right ear and blue headphone on the child’s left ear. Present a conditioning tone of 40dB at 1000 Hz.

Sweep Screening

Right: 20 dB at 1000, 2000, 4000 Hz

Switch to

Left: 20 dB at 4000, 2000, 1000 Hz

Increase to 25dB at 500 Hz

Switch to **Right:** 25 dB at 500 Hz

Present each sound 1-2 times.

Volume (dB)	25	20	20	20
Frequency (Hz)	500	1000	2000	4000
Right Ear				
Left Ear				

Minnesota Department of Health
Maternal and Child Health
Child and Teen Checkups
651-201-3650
health.childteencheckups@state.mn.us
www.health.state.mn.us

05/2024

*To obtain this information in a different format, call:
651-201-3650.*